

SereneLife
by **PYLE**



SLFTMSG25 - SLFTMSG35 - SLFTMSG45
Foot Massager
Heel, Toe & Ankle Massage Therapy

SereneLife
by **PYLE**



VISIT US ONLINE:
Have a question?
Need service or repair?
Want to leave a comment?
PyleUSA.com/ContactUs

PYLE

**PLEASE READ ALL INSTRUCTIONS AND
SAVE THESE INSTRUCTIONS FOR FUTURE REFERENCE..**

IMPORTANT SAFETY INSTRUCTIONS:


WHEN USING AN ELECTRICAL APPLIANCE, BASIC PRECAUTIONS SHOULD ALWAYS BE FOLLOWED. READ ALL INSTRUCTIONS BEFORE USING THIS APPLIANCE

- This appliance should never be left unattended when plugged in. Unplug from the outlet when not in use, and before putting on or taking off parts.
- Do not operate under a blanket or pillow. Excessive heating can occur and cause fire, electric shock, or injury to people.
- Close supervision is necessary when this appliance is used near children or disabled persons.
- Use this appliance only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
- Never operate this appliance if it has a damaged cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the appliance to a service center for examination and repair.
- Do not carry this appliance by the supply cord or use the cord as a handle.
- Keep the cord away from heated surfaces.
- Never operate the appliance with the air openings blocked. Keep the air openings free of lint, hair, and the like.
- People who are heat-sensitive should pay attention when the surface of the unit becomes hot.
- Never drop or insert any object into any opening.
- Do not use outdoors.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.

IF YOU HAVE A MEDICAL CONDITION LIKE DIABETES, SUFFER FROM BLOOD CLOTS, RECENTLY HAD SURGERY AND/OR ARE PREGNANT, CONSULT YOUR DOCTOR BEFORE USING THE UNIT

- Do not use the massager for any purpose other than massage.
- **KEEP OUT OF CHILDREN'S REACH.**
- Do not use the massager if some parts are missing.
- Do not insert any metal objects inside the massager. This may result in electric shock.
- Do not stand inside the unit or apply excessive pressure to the massager. This may significantly slow down or stop the motor, causing the unit to overheat.
- Do not fall asleep using the massager.
- Do not cover the massager during operation; this may result in overheating of the unit.
- Do not unplug the massager from the electrical outlet during operation. This may result in injury.
- Due to risk of electric shock do not place wet legs and feet inside the unit.
- Do not use the massager in the bathroom where you might come in contact with water.
- Prior to use make sure that your voltage corresponds with the voltage of the massager.
- Do not pull on the power cord when disconnecting from the electrical outlet. Hold the plug and unplug properly from outlet.
- Due to risk of electric shock, do not handle the plug and cord with wet hands.
- Do not wrap cord around the unit. This may damage the cord and result to electric shock or fire.
- Do not place power cord under the massager or any heavy object.
- Do not operate massager with a damaged cord, do not try to fix it yourself, contact a licensed electrician.
- If you experience any discomfort while using the massager, press the "Power OFF" button to turn the unit OFF.
- Some users may experience itching of feet during and after the massage. This happens as a result of improved blood circulation and is OK.


REMOTE CONTROL



USAGE:


1. Connect the provided adaptor to the massager, and then plug into the wall outlet.
2. Choose a comfortable chair of equal height. Place your legs into the slots and relax.
3. Press the Power Button to turn the massager on. The unit will beep one time and the LCD will display "15" and then will begin the countdown.
4. The massager will start operating automatically.
5. There are a total of 3 levels (P1-P3). The default setting is level 2 with Heat (Massage Mode button will remain BLUE during operation, the buttons will be flashing and the Heat will stay on).
 - If desired, press the button to move to level 3 (will flash RED) and the LCD panel will briefly display P3.
 - During the massage, feel free to change levels according to your comfort level.
To increase the level, press the button.
To lower the level, press the button.
The LCD panel will briefly display the selected level.

Removable Cover



LCD Panel

Control Panel



(1) Heat key - activates the heat function

(2) Decrease Button - press to decrease the massage intensity

(3) Increase Button - press to increase the massage intensity

(4) Speed Button- press to adjust vibration speed

(5) Massage Mode - select massage modes

(6) Power Button - turns unit ON and OFF

6. You will feel the air pressure increasing and then releasing during the duration of the massage. For optimal massage experience, you will feel the change of direction and type of massage. If you feel that a certain level is too high or too low, use the buttons to adjust the level. Refer to point 5 on the previous page. There are 3 modes massage programs, the difference of each mode as following, and default setting is mode 1

- **Mode 1:** Air pressure and rolling function work at regular motion and speed, and button with blue light.
- **Mode 2:** Air pressure works the same as mode 1 rolling function work rhythmic, and button with red light.
- **Mode 3:** Only air pressure work, rolling function stop. And button with alternate blue and red light.

7. To adjust the vibration speed, press the speed button the vibration speed can change from high to low.

8. To turn the heat off, press the heat mode button the red light will be off indicating the heat function is **OFF**.

9. To turn **OFF** the unit at any time, just press the power button. The next time you turn it back **ON**, it will go back to default setting Level 2. (P2) -15 minutes.

10. Unplug the unit when done using the massager.

CARE AND STORAGE:

- Before cleaning, always make sure to turn off the unit from control panel then unplug from the electrical outlet.
- Do not use any abrasive cleaners when cleaning the unit.
- Use a dry cloth when cleaning the massager. Never use the damp cloth to clean the control panel and massager unit.
- When not in use for long periods of time, store the massager in a dust and moisture free environment with the cord coiled. Do not wrap the cord around unit.
- Do not leave the unit under direct sunlight for extended periods of time.
- Do not store the massager unit near the heat source.
- If the massager unit malfunctions or brakes down, never try to repair it yourself, consult a licensed electrician.

SLFTMSG25	Penetrating Feet Roller Mode
SLFTMSG35	Shiatsu Pressure Therapy Mode
SLFTMSG45	Vibration Pressure Therapy Mode

FEATURES:

- Full Foot Therapy from Heel to Toe to Ankle
- Simply Slide in Your Feet for Soothing Comfort
- Air Pressure Deep Kneading for Tension Relief
- Increase Blood Flow & Circulation
- Relaxes & Relieves Muscle Pain
- Button-Activated Heat Function
- Easy Plug-in Design with Hassle-Free Operation
- Adjustable Massage Intensity Levels
- Selectable Speed Settings
- Digital Touch Buttons with LCD Display
- Safety & Power Saving Auto Turn-Off
- Removable & Washable Foot Covers
- Height Adjustable Base

WHAT'S IN THE BOX:

- Foot Massager
- (2) Removable Slip-in Foot Covers
- Remote Control
- Power Cable

TECHNICAL SPECS:

- High-Powered Massage Motors: 60 Watt
- Automatic Turn-Off: 15 minutes
- Power Cable Length: 7.5' ft.
- Power: 110-240V
- Dimensions (L x W x H):17.7" x 15.7" x 11.0" -inches
- Sold as: 1